


enVoke™ Intermediate Catheter




Product Information

eNVo™ Diameter Group	REF Catalogue Number	Description	“C” Inner Diameter	“D” Usable Length	Stiffness		
eNVo™ 44	FG-005-4011	eNVo™ Intermediate Catheter, 044, 90 cm, Soft	.044 inch (1.12 mm)	90 cm	Soft		
	FG-005-4012	eNVo™ Intermediate Catheter, 044, 90 cm, Regular			Regular		
	FG-005-4013	eNVo™ Intermediate Catheter, 044, 90 cm, Support			Support		
	FG-005-4111	eNVo™ Intermediate Catheter, 044, 105 cm, Soft		105 cm	Soft		
	FG-005-4112	eNVo™ Intermediate Catheter, 044, 105 cm, Regular			Regular		
	FG-005-4113	eNVo™ Intermediate Catheter, 044, 105 cm, Support			Support		
	FG-005-4211	eNVo™ Intermediate Catheter, 044, 115 cm, Soft		115 cm	Soft		
	FG-005-4212	eNVo™ Intermediate Catheter, 044, 115 cm, Regular			Regular		
	FG-005-4213	eNVo™ Intermediate Catheter, 044, 115 cm, Support			Support		
	FG-005-4311	eNVo™ Intermediate Catheter, 044, 132 cm, Soft		132 cm	Soft		
	FG-005-4312	eNVo™ Intermediate Catheter, 044, 132 cm, Regular			Regular		
	FG-005-4313	eNVo™ Intermediate Catheter, 044, 132 cm, Support			Support		
	FG-005-4511	eNVo™ Intermediate Catheter, 044, 150 cm, Soft		150 cm	Soft		
	FG-005-4512	eNVo™ Intermediate Catheter, 044, 150 cm, Regular			Regular		
FG-005-4513	eNVo™ Intermediate Catheter, 044, 150 cm, Support	Support					
eNVo™ 60	FG-005-5011	eNVo™ Intermediate Catheter, 060, 90 cm, Soft	.060 inch (1.52 mm)	90 cm	Soft		
	FG-005-5012	eNVo™ Intermediate Catheter, 060, 90 cm, Regular			Regular		
	FG-005-5013	eNVo™ Intermediate Catheter, 060, 90 cm, Support			Support		
	FG-005-5111	eNVo™ Intermediate Catheter, 060, 105 cm, Soft		105 cm	Soft		
	FG-005-5112	eNVo™ Intermediate Catheter, 060, 105 cm, Regular			Regular		
	FG-005-5113	eNVo™ Intermediate Catheter, 060, 105 cm, Support			Support		
	FG-005-5211	eNVo™ Intermediate Catheter, 060, 115 cm, Soft		115 cm	Soft		
	FG-005-5212	eNVo™ Intermediate Catheter, 060, 115 cm, Regular			Regular		
	FG-005-5213	eNVo™ Intermediate Catheter, 060, 115 cm, Support			Support		
	FG-005-5311	eNVo™ Intermediate Catheter, 060, 132 cm, Soft		132 cm	Soft		
	FG-005-5312	eNVo™ Intermediate Catheter, 060, 132 cm, Regular			Regular		
	FG-005-5313	eNVo™ Intermediate Catheter, 060, 132 cm, Support			Support		
	eNVo™ 72	FG-005-6011		eNVo™ Intermediate Catheter, 072, 90 cm, Soft	.072 inch (1.83 mm)	90 cm	Soft
		FG-005-6012		eNVo™ Intermediate Catheter, 072, 90 cm, Regular			Regular
FG-005-6013		eNVo™ Intermediate Catheter, 072, 90 cm, Support	Support				
FG-005-6111		eNVo™ Intermediate Catheter, 072, 105 cm, Soft	105 cm	Soft			
FG-005-6112		eNVo™ Intermediate Catheter, 072, 105 cm, Regular		Regular			
FG-005-6113		eNVo™ Intermediate Catheter, 072, 105 cm, Support		Support			
FG-005-6211		eNVo™ Intermediate Catheter, 072, 115 cm, Soft	115 cm	Soft			
FG-005-6212		eNVo™ Intermediate Catheter, 072, 115 cm, Regular		Regular			
FG-005-6213		eNVo™ Intermediate Catheter, 072, 115 cm, Support		Support			
FG-005-6311		eNVo™ Intermediate Catheter, 072, 132 cm, Soft	132 cm	Soft			
FG-005-6312		eNVo™ Intermediate Catheter, 072, 132 cm, Regular		Regular			
FG-005-6313		eNVo™ Intermediate Catheter, 072, 132 cm, Support		Support			

Dimensions

eNVo™ Diameter Group	“C” Inner Diameter	“A” Distal Outer Diameter	“B” Proximal Outer Diameter
eNVo™ 44	.044 inch (1.12 mm)	.054 inch (1.38 mm)	.056 inch (1.43 mm)
eNVo™ 60	.060 inch (1.52 mm)	.073 inch (1.86 mm)	.078 inch (1.99 mm)
eNVo™ 72	.072 inch (1.83 mm)	.082 inch (2.09 mm)	.084 inch (2.14 mm)



INTENDED USE / INDICATIONS FOR USE

The eNVoKe™ catheters are intended for the introduction of interventional devices and infusion of diagnostic or therapeutic agents into the neuro, peripheral, and coronary vasculatures.

DEVICE DESCRIPTION

The eNVoKe™ catheters are variable stiffness, single lumen catheters designed to access small, tortuous vasculature. They are available in a variety of lengths, stiffnesses, and inner and outer diameters. The outer surface of the catheter has a hydrophilic coating on the distal end to reduce friction in the vessel. The catheter also incorporates a liner to facilitate introduction of interventional devices through its lumen. The proximal end has a Luer hub for attachment to other devices and the distal tip has radiopaque marker(s) to aid visualization and positioning under fluoroscopy.

COMPATIBILITY

Catheter	Minimum Guide Catheter Inner Diameter	Maximum Guidewire Diameter
eNVoKe™ 44	.088 inch (2.24 mm)	.035 inch (0.89 mm)
eNVoKe™ 60	.088 inch (2.24 mm)	.035 inch (0.89 mm)
eNVoKe™ 72	.088 inch (2.24 mm)	.035 inch (0.89 mm)

- Use of the eNVoKe™ in a lumen less than .003 inch (0.08 mm) larger than the outer diameter of the eNVoKe™ is not recommended.
- Limited testing has been performed with solutions such as contrast media and saline. Delivery of solutions other than the types tested is not recommended.
 - Compatibility with glue or glue mixtures has not been established.
 - Compatibility with DMSO or agents suspended in DMSO has not been established.
 - Ensure embolic material compatibility with eNVoKe™ prior to use.

CONTRAINDICATIONS

- None known.

WARNINGS

- Do not use if damage to the device is observed.
- Do not use if the product sterile barrier system or its packaging is compromised.
- This device is intended for single use only. Do not resterilize and/or reuse in multiple patients. Structural integrity, sterility and/or function may be impaired by resterilization or re-use.
- Never advance or withdraw an intravascular device against resistance until the cause of resistance is determined by fluoroscopy. Movement of the device against resistance could dislodge a clot, perforate a vessel wall, or damage the device.
- The catheter should be manipulated under fluoroscopy only. Do not attempt to move the catheter without observing the resultant tip response.
- Do not exceed 2070 kPa (300 psi) maximum recommended infusion pressure. Excess pressure may result in catheter damage or patient injury. Use of power injectors requires careful monitoring of catheter tip placement in the vasculature to avoid vessel damage.
- Do not exceed 317 kPa (46 psi) if lumen of catheter is occluded as this may result in catheter damage or patient injury.

POTENTIAL COMPLICATIONS

Possible complications of the use of the eNVoKe™ catheter include, but are not limited to:

- Adverse reaction to antiplatelet/anticoagulation agents or contrast media
- Additional surgical intervention
- Allergic reaction
- Aneurysm or pseudoaneurysm
- Aneurysm perforation or rupture
- Anesthesia complications
- Arteriovenous fistula
- Cerebral infarct
- Change in mental status
- Death
- Device deformation, collapse, fracture, or malfunction
- Edema, including braid and pulmonary
- Embolization (air, tissue, or thrombotic emboli)
- Hematoma at the site of entry
- Hemorrhage
- Hypotension/hypertension
- Infection
- Inflammation, including sterile inflammation or granulomas at the access site
- Intracranial hemorrhage including subarachnoid hemorrhage and hemorrhagic transformation
- Ischemia
- Myocardial embolism
- Myocardial infarct
- Neurological deficits including stroke and death
- New stroke/ cerebrovascular accident/ transient ischemic attack (TIA)
- Pain at the site of entry
- Post procedural bleeding
- Pulmonary embolism
- Pulmonary infarct
- Pseudoaneurysm
- Renal failure
- Respiratory failure
- Shock
- Stroke
- Thrombosis (acute and subacute)
- Tissue necrosis, transient or long-lasting
- Vascular thrombosis
- Vascular occlusion
- Vasospasm (Vasospasm)
- Vessel trauma, dissection, perforation, rupture, or injury

eNVoKe™ COMPATIBILITY

Non-clinical compatibility testing has been performed with the eNVi™-SR Retriever.^{1,2}

WARNINGS (CONTINUED)

- The safety and effectiveness of the device has not been established, or is unknown, in vascular regions other than those specifically indicated.
- This device is coated with a hydrophilic coating at the distal end of the device for a length indicated on the label. Please refer to Procedure section for further information on how to prepare and use this device to ensure it performs as intended. Failure to abide by the warnings in this labeling may result in damage to the device coating, which may necessitate intervention or result in serious adverse events.

PRECAUTIONS

- Carefully read these directions before using this product. Observe warning and safety precautions.
- The appropriate anti-platelet and anti-coagulation therapy should be administered in accordance with standard medical practice.
- The device should only be used by physicians experienced in angiographic and percutaneous interventional procedures, at medical facilities with the appropriate fluoroscopic equipment.
- Use device prior to "Use-by" date printed on label.
- Prevent exposure to temperatures in excess of 55°C. Exposure to temperatures above this temperature may damage device and accessories. Do not autoclave.
- To prevent thrombus formation and contrast media crystal formation, maintain a constant infusion of appropriate flush solution between guide catheter and microcatheter and between the microcatheter and guidewire.
- If flow through catheter becomes restricted, do not attempt to clear catheter lumen by infusion. Doing so may cause catheter damage or patient injury. Remove and replace catheter.
- Torquing the catheter may cause damage which could result in kinking to separation of the catheter shaft.
- If an intraluminal device becomes lodged in the catheter, or if the catheter becomes severely kinked, withdraw the entire system (intraluminal device, catheter, and guiding catheter).
- The Shaping Mandrel is not intended for use inside the body. Ensure Shaping Mandrel is removed from the catheter prior to introduction into the RHV or other accessories.
- Use only a steam source to shape the catheter tip. Do not use other heat sources or the catheter may be damaged.
- Verify that the diameter of any guidewire or accessory device that is used is compatible with the inner diameter of the catheter prior to use.
- The eNVoKe™ catheter has a lubricious hydrophilic coating on the outside of the distal end.
 - The coating must be kept hydrated in order to be lubricious.
 - The coating is incompatible with solvents such as alcohols or cleaning agents. Avoid using alcohols, antiseptic solutions, or other solvents as these may damage the coating, which could affect device safety and performance.
 - Avoid wiping the device with dry gauze as this may damage the device coating. Avoid excessive wiping of the coating.
 - Avoid device insertion through a metal cannula or needle. Manipulation, advancement, and/or withdrawal through a metal device may result in coating material remaining in the metal device leading to adverse events such as embolization.
- Caution: Federal (USA) law restricts this device to sale, distribution and use by or on the order of a physician.



EC REP
Emergo Europe
 Prinsessegracht 20
 2514 AP The Hague
 The Netherlands

CE 0297
 Pending FDA
 Clearance, Not
 Available for Sale in
 the United States.

NeuroVasc Technologies, Inc.
 3 Jenner, Suite 100
 Irvine, CA 92618, USA
 +1.949.258.9946

¹TR20-001, ²TR20-003
PROM-018 Rev. A (2022-04)



NeuroVascTechnologies.com