

WEGO

LEADING BRAND OF MEDICAL
DEVICES IN CHINA



WEIGAO MEDICAL INTERNATIONAL CO., LTD.

ADD: No.1 Weigao Road,Hi-Tech Zone,Weihai,Shandong,China
TEL: +86 631 5660699 FAX: +86 631 5622563 P.C.: 264210
E-mail: info@weigaogroup.com exp@weigaogroup.com
Http: //www.wegomedical.com/

This copyright belongs to Wego Medical International Co., Ltd. Our company reserves the right to investigate legal liability in case of infringement.

2024.2 EDITION





Disposable Surgical Intervention
Auxiliary Package

WG04-006-001

For surgical use. Components contain surgical gowns, sheets, drapes, fenestrated drapes and disposable plastics.



Surgical Package for Single Use

WG04-006-002

For surgical use. Suitable for various surgeries like orthopedic, endoscopy, ophthalmic, neurosurgery. Components contain surgical gowns, sheets, drapes, fenestrated drapes, and disposable plastics.



Vein Puncture Assistant Kit

WG04-006-003

Combination of dressing and adhesive film



Puncture Kit for Single Use

WG04-006-004

Combined needle eases the procedure of puncture and drainage



Blood Transfusion
Consumables



Single bag is used for collection of whole blood with anticoagulant CPDA-I solutions USP.

This Single blood bag system with a primary bag contains anticoagulant CPDA-I (citrate-phosphat-Dextrose-Adenine) which preserves red blood cells up to 35 days at 2°C-6°C.

Sizes and Specification(WG02-001-001)

Code Number	Bag Capacity	Anticoagulant	Needle	Al foil pack /carton	
SC-100W	100ml	CPDA-I	14ml	16G	5/100
SC-250W	250ml	CPDA-I	35ml	16G	5/100
SC-350W	350ml	CPDA-I	49ml	16G	5/100
SC-450W	450ml	CPDA-I	63ml	16G	5/100
SC-500W	500ml	CPDA-I	70ml	16G	5/100

Code Number	Bag Capacity	Anticoagulant	Needle	Al foil pack /carton	
SC-100E	100ml	CPDA-I	14ml	16G	5/100
SC-250E	250ml	CPDA-I	35ml	16G	5/100
SC-350E	350ml	CPDA-I	49ml	16G	5/100
SC-450E	450ml	CPDA-I	63ml	16G	5/100
SC-500E	500ml	CPDA-I	70ml	16G	5/100



This triple blood bag system is used for the separation of three components from whole blood. The triple system includes one primary bag with anticoagulant CPDA-I solution BP and two empty satellite bags.

The triple SAGM blood bag system includes one primary bag with anticoagulant CPD solution USP. SAGM solution in another bag and one empty satellite bag.

Sizes and Specification(WG02-001-003)

Code Number	Bag Capacity	Anticoagulant	Needle	Al foil pack /carton	Satellite	
TC-250W	250ml	CPDA-I	35ml	16G	2/50	2X200ml
TC-350W	350ml	CPDA-I	49ml	16G	2/50	2X300ml
TC-450W	450ml	CPDA-I	63ml	16G	2/50	2X300ml
TC-500W	500ml	CPDA-I	70ml	16G	2/50	2X300ml
TS - 350W	350ml	CPD+SAGM	49ml+78ml	16G	2/50	2X300ml
TS - 450W	450ml	CPD+SAGM	63ml+100ml	16G	2/50	2X300ml
TS - 500W	500ml	CPD+SAGM	70ml+111ml	16G	2/50	2X300ml

Code Number	Bag Capacity	Anticoagulant	Needle	Al foil pack /carton	Satellite	
TC-250E	250ml	CPDA-I	35ml	16G	2/50	2X200ml
TC-350E	350ml	CPDA-I	49ml	16G	2/50	2X300ml
TC-450E	450ml	CPDA-I	63ml	16G	2/50	2X300ml
TC-500E	500ml	CPDA-I	70ml	16G	2/50	2X300ml
TS - 350E	350ml	CPD+SAGM	49ml+78ml	16G	2/50	2X300ml
TS - 450E	450ml	CPD+SAGM	63ml+100ml	16G	2/50	2X300ml
TS - 500E	500ml	CPD+SAGM	70ml+111ml	16G	2/50	2X300ml



This double blood bag system is used for separation of two components from whole blood. This double system includes one primary bag with anticoagulant CPDA-I solutions USP and one empty satellite bag.

Sizes and Specification(WG02-001-002)

Code Number	Bag Capacity	Anticoagulant	Needle	Al foil pack /carton	Satellite	
DC-250W	250ml	CPDA-I	35ml	16G	5/80	1X200ml
DC-350W	350ml	CPDA-I	49ml	16G	5/80	1X300ml
DC-450W	450ml	CPDA-I	63ml	16G	5/80	1X300ml
DC-500W	500ml	CPDA-I	70ml	16G	5/80	1X300ml

Code Number	Bag Capacity	Anticoagulant	Needle	Al foil pack /carton	Satellite	
DC-250E	250ml	CPDA-I	35ml	16G	5/80	1X200ml
DC-350E	350ml	CPDA-I	49ml	16G	5/80	1X300ml
DC-450E	450ml	CPDA-I	63ml	16G	5/80	1X300ml
DC-500E	500ml	CPDA-I	70ml	16G	5/80	1X300ml



This quadruple blood bag system is used for the separation of four components from whole human blood. The quadruple system includes one primary bag with anticoagulant CPDA-I solution USP and three empty satellite bags.

The quadruple SAGM blood bag system includes one primary bag with anticoagulant CPD solution USP, SAGM solution in another bag and two empty satellite bags.

Sizes and Specification(WG02-001-004)

Code Number	Bag Capacity	Anticoagulant	Needle	Al foil pack /carton	Satellite	
QC-350W	350ml	CPDA-I	49ml	16G	2/40	3X300ml
QC-450W	450ml	CPDA-I	63ml	16G	2/40	3X300ml
QC-500W	500ml	CPDA-I	63ml	16G	2/40	3X300ml
QS-350W	350ml	CPDA-I	70ml	16G	2/40	3X300ml
QS-450W	450ml	CPD+SAGM	49ml+78ml	16G	2/40	3X300ml
QS-500W	500ml	CPD+SAGM	63ml+100ml	16G	2/40	3X300ml

Code Number	Bag Capacity	Anticoagulant	Needle	Al foil pack /carton	Satellite	
QC-350E	350ml	CPDA-I	49ml	16G	2/40	3X300ml
QC-450E	450ml	CPDA-I	63ml	16G	2/40	3X300ml
QC-500E	500ml	CPDA-I	63ml	16G	2/40	3X300ml
QS-350E	350ml	CPDA-I	70ml	16G	2/40	3X300ml
QS-450E	450ml	CPD+SAGM	49ml+78ml	16G	2/40	3X300ml
QS-500E	500ml	CPD+SAGM	63ml+100ml	16G	2/40	3X300ml



safety needle protector and blood sampling collection device.