


**WEGO**



CE 0123 ISO

## High-pressure Angiographic Syringe for Single Use and Accessory



**WEIGAO MEDICAL INTERNATIONAL CO., LTD.**

ADD: No.1 Weigao Road, Hi-Tech Zone, Weihai, Shandong, China  
TEL: +86 631 5660699 FAX: +86 631 5622563 P.C.: 264210  
E-mail: [info@weigaogroup.com](mailto:info@weigaogroup.com) [exp@weigaogroup.com](mailto:exp@weigaogroup.com)  
<https://www.wegomedical.com/>

This copyright belongs to Wego Medical International Co., Ltd. Our company reserves the right to investigate legal liability in case of infringement.

2024.2 EDITION



# CT SINGLE

Computed Tomography



## CT-200A-MA

**Compatible with:** Medrad Stellant  
**Highest pressure limit:** 600psi  
**Accessories:**  
CT Connecting Tube/Quick Fill Tube/Spike

## CT-200-LF

**Compatible with:** LF CT9000m, Optiventage  
**Highest pressure limit:** 600psi  
**Accessories:**  
CT Connecting Tube/Quick Fill Tube/Spike

## CT-200-NE

**Compatible with:** Nemoto-A25,A60  
Medtron Accutron  
**Highest pressure limit:** 300psi  
**Accessories:**  
CT Connecting Tube/Quick Fill Tube/Spike

## CT-200-MD

**Compatible with:** Medtron Accutron CT  
**Highest pressure limit:** 300psi  
**Accessories:**  
CT Connecting Tube/Quick Fill Tube/Spike



CT-200D-MA

**Compatible with:**Zenith-C20 B  
**Highest pressure limit:**300psi  
**Accessories:**  
CT Connecting Tube/Quick Fill Tube/Spike



CT-200-MA

**Compatible with:**Medrad MCT Plus,  
Vistron CT,Envision CT  
**Highest pressure limit:**300psi  
**Accessories:**  
CT Connecting Tube/Quick Fill Tube/Spike



CT-120-AK

**Compatible with:**ASA-200H  
**Highest pressure limit:** 145psi  
**Accessories:**  
CT Connecting Tube/Quick Fill Tube/Spike



CT-200-BI

**Compatible with:**Imaxeon Salient  
**Highest pressure limit:**300psi  
**Accessories:**  
CT Connecting Tube/Quick Fill Tube/Spike



## CT DOUBLE

Computed Tomography




CT-200+200-EZ

**Compatible with:** EZEM Empower  
**Highest pressure limit:** 600psi  
**Accessories:**  
 Y-Connecting Tube/Quick Fill Tube/Spiket



CT-200+200-NE

**Compatible with:** Nemoto Dual Shot  
**Highest pressure limit:** 300psi  
**Accessories:**  
 Y-Connecting Tube/Quick Fill Tube/Spike



CT-200+200-BI

**Compatible with:** Imaxeon Salient  
**Highest pressure limit:** 300psi  
**Accessories:**  
 Y-Connecting Tube/Quick Fill Tube/Spike




CT-50+200-NE

**Compatible with:** Nemoto Dual Shot  
**Highest pressure limit:** 300psi  
**Accessories:**  
 Y-Connecting Tube/Quick Fill Tube/Spike



CT-200+200-MA

**Compatible with:** Medrad Stellant  
**Highest pressure limit:** 600psi  
**Accessories:**  
 Y-Connecting Tube/Quick Fill Tube/Spike



CT-200+200-MD

**Compatible with:** Medtron Accutron  
**Highest pressure limit:** 300psi  
**Accessories:**  
 Y-Connecting Tube/Quick Fill Tube/Spike



CT-200+200-LF

**Compatible with:** LF Optivantage  
**Highest pressure limit:** 600psi  
**Accessories:**  
 Y-Connecting Tube/Quick Fill Tube/Spike



# MRI

Magnetic Resonance Imaging



## MRI-60+110-MA

**Compatible with:** Medrad Spectris  
**Highest pressure limit:** 200psi  
**Accessories:** Y-Connecting Tube/Quick Fill Tube/Spike

## MRI-60+60-LF

**Compatible with:** LF Optistar  
**Highest pressure limit:** 145psi  
**Accessories:** Y-Connecting Tube/Quick Fill Tube/Spike

## MRI-50+50-NE

**Compatible with:** Nemoto Sonic Shot 50:  
**Highest pressure limit:** 300psi  
**Accessories:** Y-Connecting Tube/Quick Fill Tube/Spike

## MRI-60+60-MA

**Compatible with:** Medrad Spectris  
**Highest pressure limit:** 200psi  
**Accessories:** Y-Connecting Tube/Quick Fill Tube/Spike



# DSA

Magnetic Resonance Imaging



## DSA-I50B-LF

**Compatible with:** LF Angiomat Illumena  
**Highest pressure limit:** 1200psi  
**Accessories:**  
CT Connecting Tube/Quick Fill Tube/Spike

## DSA-I50A-LF

**Compatible with:** LF Angiomat 6000  
**Highest pressure limit:** 1200psi  
**Accessories:**  
CT Connecting Tubing/Quick Fill Tube/Spike

## DSA-I50-MA

**Compatible with:** Medrad Mark V  
**Highest pressure limit:** 1200psi  
**Accessories:**  
CT Connecting Tube/Quick Fill Tube/Spike

## DSA-90-NE

**Compatible with:** Nemoto M800  
**Highest pressure limit:** 1200psi  
**Accessories:**  
CT Connecting Tube/Quick Fill Tube/Spike

# ACCESSORY



STERILE PRESSURE CONNECTING TUBE FOR SINGLE USE

Type: Y type coiled  
Valve: With check valve  
Quantity: 10pcs/box  
100pcs/carton  
Material: PVC, PP, PC  
Pressure: 600psi

STERILE PRESSURE CONNECTING TUBE FOR SINGLE USE

Type: I type coiled  
Valve: Without check valve  
Quantity: 10pcs/box  
100pcs/carton  
Material: PVC, PP, PC  
Pressure: 600psi

STERILE PRESSURE CONNECTING TUBE FOR SINGLE USE

Type: I type coiled  
Valve: With check valve  
Quantity: 10pcs/box  
100pcs/carton  
Material: PVC, PP  
Pressure: 330psi

STERILE TRANSFER SET FOR SINGLE USE

Type: Non-coiled  
Valve: Dual check valve  
Quantity: 10pcs/box  
100pcs/carton  
Material: PVC, PP, PC,ABS;Pressure: 330psi

TYPE	MACHINE MODEL	HIGHEST PRESSURE LIMIT	ACCESSORY
CT-200-MA	Medrad MCT Plus Vistron CT,Envision CT	300psi	CT Connecting Tube/ Quick Fill Tube/Spike
CT-200A-MA	Medrad Stellant	600psi	
CT-200D-MA	Zenith-C20 B	300psi	
CT-100-NE	Nemoto-A25,A60,Dual shot	300psi	
CT-200-NE	Nemoto-A25,A60	300psi	
CT-200-LF	LF CT9000, Optivantage	600psi	
CT-150-BI	Imaxeon Visimax	145psi	
CT-140-BI	Biotel PJ3 MK2	145psi	
CT-200-BI	Imaxeon Salient	300psi	
CT-120-AK	ASA-200H	145psi	
CT-200-AK	ASA-200E	145psi	
CT-200-EZ	EZEM Empower	600psi	
CT-200-MD	Medtron Accutron	300psi	

TYPE	MACHINE MODEL	HIGHEST PRESSURE LIMIT	ACCESSORY
DSA-90-NE	Nemoto M800	1200psi	CT Connecting Tube/ Quick Fill Tube/Spike
DSA-130-MA	Medrad Mark IV	1200psi	
DSA-150-MA	Medrad Mark V	1200psi	
DSA-120-NE	Nemoto 120S	1200psi	
DSA-150A-LF	LF Angiomat 6000	1200psi	
DSA-150B-LF	LF Angiomat Illumena	1200psi	

TYPE	MACHINE MODEL	HIGHEST PRESSURE LIMIT	ACCESSORY
CT-200+200-MA	Medral Stellant	600psi	Y-Connecting Tube/ Quick Fill Tube/Spike
CT-200+200-MD	Medtron Accutron	300psi	
CT-50+100-NE	Nemoto Dual Shot	300psi	
CT-50+200-NE	Nemoto Dual Shot	300psi	
CT-200+200-NE	Nemoto Dual Shot	300psi	
CT-100+100-NE	Nemoto Dual Shot	300psi	
CT-200+200-LF	LF CT9000,Optivantage	600psi	
CT-200+200-BI	Imaxeon Salient	300psi	
CT-200+200-AK	ASA-300	145psi	
CT-200+200-EZ	EZEM Empower	600psi	

TYPE	MACHINE MODEL	HIGHEST PRESSURE LIMIT	ACCESSORY
MRI-100+100-EZ	EZEM EmpowerMR	300psi	Y-Connecting Tube/ Quick Fill Tube/Spike
MRI-60+110-MA	Medrad Solaris	200psi	
MRI-60+60-MA	Medrad Spectris	200psi	
MRI-60+60-LF	LF Optistar	145psi	
MRI-60+60-MD	Medtron Accutron MR	300psi	
MRI-50+50-NE	Nemoto Sonic Shot 50	300psi	